

# Protocol for Zymogram

2-17-06

Lee, Sunyoung

## Reagents:

1% Gelatin in H<sub>2</sub>O (Fisher Blood 275)

1% Casein (Sigma)

SDS-PAGE gel stock w/o urea

Wash buffer: 2.5% Triton X-100 in H<sub>2</sub>O (+0.02% NaN<sub>3</sub>)

Incubation Buffer: 50mM Tris-HCl (pH8.0), 5mM CaCl<sub>2</sub>, 0.02% NaN<sub>3</sub>

## Gels:

Regular separating gel containing 10-12% substrate

Regular stacking gel

## Protocol:

1- Collect media from cells

(if desired, inactivate non-metalloproteases with PMSF and/or NEM)

2- Centrifuge to remove cellular debris

(if necessary concentrate with centricon units or dialyze and lyophilize)

3- Add Laemmli loading buffer (**OMIT UREA AND REDUCING AGENTS, DO NOT HEAT**)

4- Load samples directly onto gel

5- Run gel

6- Wash 2X 20min in wash buffer

7- Wash 10min in incubation buffer

8- Place gel in sealable container with fresh incubation buffer and incubate at 37°C for 24h to 48h

9- Fix and stain with fresh Coomassie Blue solution

10-Destain with MeOH:AcOH:H<sub>2</sub>O(5:1:5)

11-Replace with 10% AcOH and continue destaining

12-Photograph and dry gel for storage