Immunohistochemistry Protocol for Whole-Mount or Vibratome-Sectioned Tissue Notch Antibody Staining

Shahla Jilani

For use with mouse tissue (whole-mount or vibratome-sectioned), fixed with PFA (2-4%)

- 1. Block with 0.3% Triton/1X PBS solution containing 5% Normal Donkey Serum (Jackson Immunoresearch, Catalog # 017-000-121) [or whichever species your secondary antibody is made in]
 - -2 hours
 - -room temperature
- 2. Incubate with rabbit anti-human Cleaved Notch 1 (Val1744) antibody (Cell Signaling Technology), diluted 1:250-1:500 (depending on tissue) in blocking buffer/ OR rabbit anti-rat 5261-P Notch 1 ECD antibody diluted 1:1000 in blocking buffer
 - -12-15 hours (overnight)
 - $-4^{\circ}C$
- 3. Wash with 1X PBS
 - -1 hour
 - -room temperature
- 4. Incubate with Cy3 (or FITC) –labeled donkey secondary anti-rabbit IgG antibody (Jackson Immunoresearch, Catalog # 711-165-152), diluted 1:200 in blocking buffer
 - -3-4 hours
 - -4^{0} C
- 5. Wash with 1X PBS
 - -2 hours-overnight (best if overnight)
 - $-4^{\circ}C$