

## Transfection of [T47D] Cells

### Day 1 (8/20/02):

- 1) Replace media with fresh media (10% FBS + Pen/St).
- 2) Split 1 plate's worth into a 6-well plate  
e.g.- 3mls cell suspension into each 10-cm plate  
-500 $\mu$ l cell suspension into each well of the 6-well plate
- 3) Incubate for 3-4 hours at 37°C
- 4) Dilute the adenovirus as follows:

Tube #	Virus	DMEM*	Dilution
1	5ul stock	45ul	1:10
2	5ul Tube 1	5ml	1:10 <sup>5</sup>
3	500ul Tube 2	5ml	1:10 <sup>6</sup>
4	500ul Tube 3	5ml	1:10 <sup>7</sup>
5	500ul Tube 4	5ml	1:10 <sup>8</sup>
6	500ul Tube 5	5ml	1:10 <sup>9</sup>

\*DMEM F12 (50/50 Mix) + 10%FBS + Pen/St

- 5) Add 1.5 mls of each virus/DMEM solution to each well as follows:

<b>Control (no virus)</b>	<b>1:10<sup>5</sup></b> (Tube 2)	<b>1:10<sup>6</sup></b> (Tube 3)
<b>1:10<sup>7</sup></b> (Tube 4)	<b>1:10<sup>8</sup></b> (Tube 5)	<b>1:10<sup>9</sup></b> (Tube 6)

- 6) Incubate overnight.

### Day 2 (8/21/02):

- 1) Remove media
- 2) Add SERUM-FREE DMEM and incubate for 30 minutes.
- 3) Prepare a solution of serum-free DMEM containing 5 $\mu$ g/ml heparin  
--1.5ml will be needed for each well, so make ~10ml
- 4) Remove the media
- 5) Add 1.5ml DMEM + heparin solution to each well
- 6) Incubate overnight.

### Day 3 (8/22/02):

Harvest the media.

